

Filtear® Mono

STERILE
LUBRICANT
EYE DROPS

Macroglol 400 (Polyethylene Glycol 400) BP
0.4% & Propylene Glycol BP 0.3%

Presentation

Filtear® Mono Sterile Lubricant Eye Drops: Each single use container contains Macroglol 400 (Polyethylene Glycol 400) BP 4 mg/ml & Propylene Glycol BP 3 mg/ml.

Description


Filtear® Mono contains Macroglol 400 (Polyethylene Glycol 400) (0.4%) & Propylene Glycol (0.3%). Polyethylene Glycol 400 (0.4%) & Propylene Glycol (0.3%) are eye lubricants. Eye lubricants keep the eye moist, helps to protect the eye from injury and infection, and decrease symptoms of dry eye such as burning, irritation and feeling as if something is in the eye. It is scientifically formulated to shield eyes from dry eye discomfort so that eyes feel moist and refreshed longer.

Indications and uses

It is indicated for the temporary relief of burning and irritation due to dry eye.

Dosage and Administration

To open the single use container, first hold the unit by the label end, remove a single dose unit from the pack (see picture 1). Twist and pull the head of the container to remove (see picture 2). Instill 1 or 2 drops in the affected eye(s) as needed or as directed by the physician (see picture 3). While closing the container, turn the head at a 90° angle in relation to the container, and with a slight force, press the head onto the container tip. (see picture 4). After each use close the head properly. Discard the single use container after 12 hours of opening.



Side effects

Generally well tolerated. A very serious allergic reaction to this drug is rare. It should not be used if allergic condition occurs to any ingredients of the product.

Precautions

Concomitant ocular medication should be administered 15 minutes apart to the instillation of this eye drop. To avoid contamination, do not touch tip of container to any surface. Do not use after 12 hours of opening. Do not touch single use container tip to eye. If the solution changes color or becomes cloudy, do not use the eye drops.

Contraindication

It is contraindicated in patients with known sensitivity to any ingredient in this formulation.

Use in pregnancy and lactation

Use in Pregnancy: There is no data on the use of Polyethylene Glycol 400 (0.4%) & Propylene Glycol (0.3%) during pregnancy and lactation in human. Animal studies did not show harmful effects with the active ingredient Polyethylene Glycol 400 (0.4%) & Propylene Glycol (0.3%).

Use in Lactation: Polyethylene Glycol 400 (0.4%) & Propylene Glycol (0.3%) is not absorbed systemically; there is no known potential for excretion in human breast milk.

Pediatric use

The safety and effectiveness in pediatric patients have not been established.

Drug interaction

No information is available.

Overdose

No information is available.

Storage

Do not store above 30°C. Keep away from light and out of the reach of children. Discard the single use container after 12 hours of opening.

Commercial Pack

Filtear® Mono Sterile Lubricant Eye Drops: Each single use container contains 0.5 ml of Sterile Lubricant Eye Drops.

Manufactured by

 **Incepta Pharmaceuticals Ltd**

Dhamrai Unit, Dhaka, Bangladesh

® Registered Trademark